LAMPADA CHIRURGICA ODONTOIATRICA
SURGICAL DENTAL LIGHT

ORION Tech DS

ORION Tech DS

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Illuminated field:
**19x12 cm rectangular field: light intensity
from 0 to 60.000 lux at 70 cm distance.**
**27 cm circular field: light intensity from
40.000 to 60.000 lux at 100 cm distance.**

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
Why use a surgical light?

During interventions, illumination must transmit precise information about the working area to the surgeon, while ensuring maximum eye protection. Orion Tech unique technology can obtain these 2 purposes.

Patented indirect lighting.

Orion Tech integrates the patented technique of indirect lighting. 30 inverted LEDs light 30 single domes, which offer maximum scialytic effect and guarantee an incredible illumination of the back of the mouth even in the presence of instruments. This unique technique permits working in a standing position under the light, with illumination of the oral cavity.

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Orion Tech is an evolution of the well-known Orion 40 DS light, designed to meet the requirements of the implantologist.

Orion Tech is a surgical light developed specifically for dentistry, thanks to the features tailored to the dentist's needs: at 70 cm, the rectangular illuminated field of 19 x 12 cm offers a light intensity from 0 to 60.000 lux with two colour temperatures of 4.500°K or 5.000°K. Orion Tech is a multifunctional light which can vary the size of the illuminated field, giving the possibility of switching from the dental field (19x12 cm) to the surgical field (27 cm, circular) with a simple click. Even the surgical field permits operation at double colour temperature of 4.500°K or 5.000°K and a light intensity from 40.000 to 130.000 lux.

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Double joint anti-swinging system.

Orion Tech offers a unique positioning quality thanks to a very steady arm, which maintains its position without swinging. The unique, easy Orion Tech manoeuvrability of the single joint model is easily converted to the double joint model, offering maximum visibility for all operators involved in 4/6 hand work.

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E-Dentalight

This function lights two 6 LED modules, producing a 19x12 cm illuminated field of 60.000 lux light intensity.

Cold light

Thermal irradiation increase in the working area is proportional to the number and power of the LEDs used. Thanks to the indirect light principle, Tecno-Gaz uses the lowest LED number with cold light.

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All rights reserved. The images and contents are subject to changes without prior notice. Tecno-Gaz S.p.A. cannot be held responsible for any damages resulting from incorrect or incomplete information contained herein. The images used in this document are for illustrative purposes only.

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