

 \odot

Dr's Cam V

280° twist head

1 Features of V-series

Auto focus Can get clear image with auto focus function in oral cavity.

FHD image Can get high definition image without distortion.

AF

淡:

Built-in mouse (Dr's cam only) Useable as a mouse.

Adjustable LED brightness To avoid image quality loss by light reflection, can adjust LED brightness.

Color impressions Adjust the color impression on camera to fit different monitors or screen.

280° twist head Easy to get approach anywhere in oral cavity.

Black & White function Can take and send analog X-ray film image to PC.

Plug & Play No need additional driver to install the device.

Easy to capture Easy to image capture by camera only without extra device like foot pedal.

Software compatibility

Can be used camera with most of dental imaging software. * It would need twain driver for some of dental imaging software.

2	Minimum	System	Requirement
---	---------	--------	-------------

PC system:	3.3GHz CPU, 4GB RAM		
OS:	Windows 7		
USB port:	USB 2.0		



Ez SHOT HD V

Extra Illumination Optic System

EIOS

 \odot

NEW ARRIVAL

۲

۲

۲





Fire

DrsViewer2 Image Viewer included.



octors Germany GmbH 1-Domagk-Straße 2

53121 Bonn GERMANY

Tel.: + 49 (0)228 - 53 44 14 65 Fax: +49 (0)228-53 44 14 66 www.gooddrs.de info@gooddrs.de

۲

 \bigcirc

٢