

Dental Camera

-IC-WHT60 -IC-WHW61
 -IC-DCAM70 -IC-DCAU71
 -IC-EZAM80 -IC-EZAU81
 -IC-WHMA90 -IC-WHCD100

Curing Light

-CL-DLT30 -CL-DL10
 -CL-DC20 -CL-AT24


Apex Locator

-AL-DF10 -AL-DFA20



Zoom function

Ability to enlarge the image with zoom feature to support user convenience.



* Need to install "Zoom function" program on the PC.



Software compatibility

Able to be used with almost any dental imaging software in the market.



1 Specification

Model	IC-WHMA90	
Camera	Power source	5V (USB power)
	Size	Approx. 128 x 21.5 x 32 mm
	Weight	Approx. 61g (2.2oz)
	Sensor type	5M CMOS sensor
	Lighting	White LED x 12 pcs
	Focusing	Auto-focus
	Focus range	Macro to infinity
USB cable	cable length	Approx. 3 M (9.8ft)

2 Recommended PC system

PC system:	3.3GHz CPU, 4GB RAM
OS:	Windows 7/8/10
USB port:	USB 2.0




GoodDrs



Whicam M

Macro to Infinity
 Autofocus Dental Camera with Mouse



Good Doctors Germany GmbH
 Gerhard-Domagk-Straße 2
 53121 Bonn, Germany
 Tel: +49 (0) 228 53 44 14 65 Fax: +49 (0) 228 53 44 14 66
 E-mail : info@gooddrs.de Website : www.gooddrs.de



Nov.2018 rev.0

Features of Whicam M


Image quality

- Macro to infinity**
Camera is able to focus from macro mode (close) to infinity that includes a full arch and face shot.
- Auto-focus**
Camera is able to focus in any position not only inside of oral cavity but also full arch and face shot.
- FHD image support**
The device shows FHD (1080p) in live image supported by a 5M sensor.
- Color impressions**
Adjust the color impression on camera to fit different monitors or screen.
- Black&White function**
Capture and save X-ray film image to PC



Ergonomic design

- Dental handpiece design**
Most familiar design for dentists.
- Metal body case**
Metal body case increases product durability.
- Slim head size**
Slim head design makes it easier to approach the oral cavity with camera.
- Retraction head design**
Retract patient's tongue and cheek to get better image.
- Efficient brightness**
12 pcs LEDs makes for efficient brightness.
- Dual capture buttons**
Capture buttons are placed on the top and bottom side of camera to make for an easy capture in any position
- LED indicators**
LED indicators show the device status and mode.



Convenience

- Mouse function**
Built-in mouse function for user convenience.
- Zoom function**
Ability to enlarge the image with zoom feature to support user convenience.
- LED brightness adjust**
Ability to adjust the LED brightness.
- Support twain driver**
Support twain driver allows for use in almost any dental imaging software.
- Software compatibility**
Able to be used with almost any dental imaging software in the market.
- Detachable USB cable**
Detachable USB cable makes it easy to move the camera to other dental chairs.



Macro to Infinity

Camera is able to focus from macro mode (close) to infinity that includes a full arch and face shot.



FHD image support

The device shows FHD (1080p) in live image supported by a 5M sensor.



Full HD


Original image

HD



Color impressions

Adjust the color impression on camera to fit different monitors or screen.



Same color chart image with xrite

