

Whicam Story3 CS



-IC-WHT60 -IC-WHC50

-CL-AT24

-AL-DF10

Enjoy Wireless & Auto Focus!

Wireless Autofocus Intraoral Camera with Mouse

1 Features of Whicam Story3 CS

Wireless type only $\left(\begin{pmatrix} \begin{pmatrix} & & & \\ & & & \end{pmatrix} \end{pmatrix} \right)$

High speed digital wireless

High speed digital wireless provide more clear live image without interruption



Wireless receiver

Compact wireless receiver without any power supply.



Big capacity battery

Long time use with big capacity battery.



Auto focus

Real Autofocus function shows clear images.



High definition show more detail images.



Adjustable head angle

User can use the camera with prefer angle.



Gyro mouse built-in

Useable as a mouse.



280° twist head

Easy to get approach anywhere in oral cavity.



OLED screen

Bright and clear OLED screen indicates various device status.



Adjustable LED brightness

Can adjust LED brightness to avoid light reflecting.



3 Color Impressions

Can change live image color impression (Original, Warm, Cool)



Software compatibility

Support keyboard short-cut data changing program and twain driver.

2 Specification

	Model name	IC-WHT60(Wireless model)	
n.	Camera (included head)	Operation voltage	3.7V
		Measurement	Approx. 266 x 31 x 37 mm
		Weight	Approx. 140g
		Resolution	1.3 Mega pixel
		Sensor	CMOS
	Charger	Operation voltage	5V
		Measurement	Approx. 169 x 148 x 89 mn
		Weight	Approx. 170g
	Disposable cover	Dimensions	Approx. 32 x 100 mm
		Quantity	Approx. 100 pcs by weight x 2ea
		Raw material	Polyethylene

Model name	IC-WHW61(Wired mo	IC-WHW61(Wired model)	
Camera (included head)	Operation voltage	5V	
	Measurement	Approx. 210 x 31 x 37 mm	
	Weight	Approx. 75g	
	Resolution	1.3 Mega pixel	
	Sensor	CMOS	
USB cable	Length	Approx. 3M	
Disposable cover	Dimensions	Approx. 32 x 100 mm	
	Quantity	Approx. 100 pcs by weight x 2ea	
	Raw material	Polyethylene	



3 Minimum System Requirement

PC system:	2.0 Ghz CPU, 1GB RAM
OS:	Windows XP service pack 3
	(Wireless type doesn't support Windows XP 64 bit.)
USB port:	USB 2.0





Image Viewer included.











Main View

Patient Registration

Print View

53121 Bonn **GERMANY**

Tel.: +49 (0)228 -53 44 14 65 Fax: +49 (0)228 -53 44 14 66 www.gooddrs.de info@gooddrs.de

